

NITRILE POWDERLESS GLOVE



CONTENT

- 1- ABOUT US
- 2- PRODUCT SPECIFICATION
- 3- FDA CERTIFICATION
- 4- CE CERTIFICATION
- 5- EU DECLARATION OF CONFORMITY
- 6-TEST REPORT
- 7-ISO 13485:2016 & 9001:2015
- 8- GLOVE AND PACKAGING IMAGE

ABOUT US

PRODUCT SPECIFICATION

NITRILE GLOVES SPECS

POWDER FREE

Functional Benefits

- Protection from unwanted and dangerous substances
- Beaded cuff ensures easy donning and prevent roll down
- · Superior strength with better puncture resistance
- Full textured or Finger Tip textured enhances wet and dry grip
- · Thinner gauge improves tactile sensitivity
- Custom design enhances comfort and fit
- Provide an alternative solution for individuals who are allergic to natural rubber latex

Product Specifications

Material	Synthetic Nitrile latex.
Туре	Non-Sterile Powder-Free. Ambidextrous; Finger Tip Textured; Beaded Cuff; White or Coloured (Blue, Light Blue,)
Quality Standards	Conforms to ASTM D6319 Manufactured under ISO9001: 2008, ISO 13485:2003. ISO 22000:2005 Quality Management System. Manufactured from 100% nitrile (Acrylonitrile-Butadiene)
Gloves Size	Extra-small, Small, Medium, Large, Extra-large. Marked in the check box on the shipping carton with black ink.
Storage	Store in a dry and cool place, the temperature not higher than 38°C.
Shelf-fife	3 years from the date of manufacturing.

PHYSICAL DIMENSIONS

DIMENSIONS	Standards		
DIMENSIONS	ERBEN	ASTM D6319	
Length (mm)	230 min	220 min (XS, S) 230 min (M, L, XL)	
	75 ± 5 (XS)	70 ± 10 (XS)	
	85 ± 5 (S)	$80 \pm 10 (S)$	
Width (mm)	95 ± 5 (M)	95 ± 10 (M)	
Thickness-	$105 \pm 5 (L)$	$110 \pm 10 (L)$	
	115 ± 5 (XL)	120 ± 10 (XL)	
	Fingers: 0.08 mm min	Fingers : 0.050 mm min	
Single wall (mm)	Palm: 0.06 mm min	Palm : 0.05 mm min	

PHYSICAL PROPERTIES AND BIOCOMPATIBILITY

Tensile	Tensile strength (MPA) Before aging: 18Mpa min After aging: 20Mpa min	Tensile strength (MPA) Before aging: 14Mpa min After aging: 14Mpa min	
	Elongation at break (%) Before aging: 600% min After aging: 500% min	Elongation at break (%) Before aging: 500% min After aging: 400% min	
Powder Content	2 mg/glove maximum		
Protein Content	Free Protein		

Minimum Standard of Glove Specification:

MATERIAL	CUFF THICKNESS	PALM THICKNES S	FINGER THICKNES S	LENGT H	WIDTH meaning Circumfere nce (MM)
Nitrile: 100% Latex Free and DEHP Free	3,3 mil +/- 0.5 mil	4.25 mil +/- 0.5 mil	5.3 mil +/- 0.5 mil	9.65" +/-0.25"	177.8 (XS) 203.2 (S)
	0.084mm +/- 0.013mm	0.108mm +/- 0.013mm	0.135mm +/- 0.013mm	240mm	228.5 (M) 254 (L) 279.4 (XL) 304.8 (XXL)

- 100% Nitrile disposable gloves,
- Colour preference blue,
- 100 (One Hundred) gloves/pieces per box;
- Latex free, powder free;
- Silicone free chlorinated;
- Rolled cuff;
- EN455 standard; ISO
- CE Conformity; and
- FDA 510(k).

Nitrile Examination Gloves

Product Information

Disposable Nitrile Examination Gloves, Powder Free, Non-sterile

Primary Material : Acrylonitrile Butadiene Synthetic Rubber Latex

Latex Protein Content : Latex-free

Color : White/ Black/ Blue

Size : Extra-small, Small, Medium, Large And Extra-large
Design And Feature : Ambidextrous, Straight Fingers, Finger-textured,

Beaded Cuff, Polymer Coated/ Online Chlorinated

Packing : 100 Pieces Gloves Per Dispenser, 10 Dispensers Per Carton







Medical Use

Food Contact Use

General Use

Physical Dimensions

Dimono			Standards
Dimensions Length (mm)		Supérieur	ASTM D6319 / EN 455
		240 min	230 min / 240 min
	Extra-Small	75 ± 10	75 ± 10
	Small	80 ± 10	80 ± 10
Width (mm)	Medium	95 ± 10	95 ± 10
	Large	110 ± 10	110 ± 10
	Extra-large	≥ 110	≧ 110
Thickness- single wall (mm)	Palm	Min. 0.05	Min. 0.05
	Finger	Min. 0.05	Min. 0,05

Physical Properties

Criteria		Supérieur	ASTM D6319 / EN 455
Florestian (0/)	Before Aging	Min. 500	Min. 500
Elongation (%)	After Aging	Min. 400	Min. 400 / Min. 500
Tensile Strength	Before Aging	Min. 18 (MPa)	Min. 14 (MPa) / 6.0 N
rensne strengtn	After Aging	Min. 14 (MPa)	Min. 14 (MPa) / 6.0 N













Manufactured under ISO 9001 & ISO 13485 Quality Systems



PACKAGING SPECIFICATION

SIMPLY SUPÉRIEUR

1. NITRILE EXAM GLOVE POWDER FREE

Product Name	Nitrile Exam Gloves Powder Free	
Place of Origin	Manufatured in Vietnam	
Material	100% Nitrile (synthetic rubber)	
Application	Medical use, food contact use, general us	
	etc	
Brand Name	Simply Supérieur	
	(Ever Global V	ietnam Group)
Size	XS, S, N	И, L, XL
	(size is marke	ed on the box)
Туре	Disposable, Powder Free.	
	Use for both hands, surface of rough o	
	fingertips, car wrist trim.	
Color	Blue / Black / White	
Weight	3.0g/3.5g/4.0g/4.5g	
Standard	FDA 510K, CE, EN ISO 374-5:2016,	
	EN420:2003+A1:2009, EN374-1:2016	
Packing Detail	100pcs/box, 10 boxes/carton,	
	1000pcs/carton	
Supply Ability	5,000,000 cartons/month	
Minimum Order	100,000 cartons	
Specification for Carton	Size	360x260x240 mm
	Weight	~4kg/carton

FDA CERTIFICATION



September 6, 2019

Ever Global (Vietnam) Enterprise Corporation % Albert T.w. Li
Official Third Party Correspondent
Center For Measurement Standards Of Industrial
Bldg. 16, 321 Kuang Fu Rd,sec2
Hsinchu, Tw

Re: K190403

Trade/Device Name: Disposable Powder Free Nitrile Examination Glove, White Color, Tested for Use

With Chemotherapy Drugs, Disposable Powder Free Nitrile Examination Glove,

Blue Color, Tested For Use With Chemotherapy Drug

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-powdered patient examination glove

Regulatory Class: Class I Product Code: LZA, LZC Dated: August 27, 2019 Received: September 3, 2019

Dear Albert T.w. Li:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/efdocs/efpmn/pmn.efm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of

K190403 - Albert T.w. Li

Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Elizabeth Claverie, M.S.
Assistant Director for THT4B2
Acting Assistant Director for THT4B1
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

No. 1 Road, Long Thanh Industrial Zone, Taman Village, Long Thanh District, Dong Nai Province, Vietnam
Tel: 84-61-3514022

DEPARTMENT OF HEALTH AND HUMAN SERVICES. Form Approved: OMB No. 0910-0120 Food and Drug Administration Expiration Date: 06/30/2020 Indications for Use See PRA Statement below. 510(k) Number (if known) K190403 Device Name Disposable Powder Free Nitrile Examination Glove, Blue Color, Tested For Use With Chemotherapy Drugs Indications for Use (Describe) A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs: Test Chemotherapy Drug and Concentration Minimum Breakthrough Detection Time (Min) 1. Carmustine (BCNU), 3.3 mg/ml 62 2. Thiotepa, 10.0 mg/ml 38.8 3. Cyclophosphamide (Cytoxan), 20.0 mg/ml > 240 4. Cisplatin, 1.0 mg/ml > 2405. Doxorubicin Hydrochloride, 2.0 mg/ml > 2406. Fluorouracil, 50.0 mg/ml ≥ 240 7. Dacarbazine (DTIC), 10.0 mg/ml ≥ 240 8. Etoposide (Toposar), 20.0 mg/ml ≥ 240 9. Paclitaxel (Taxol), 6.0 mg/ml ≥ 240 The maximum testing time is 240 minutes. Please note that the following drug has an extremely low permeation time: Carmustine (BCNU), 3.3 mg/ml 6.2 minutes Thiotepa, 10.0 mg/ml 38.8 minutes Type of Use (Select one or both, as applicable) Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

FORM FDA 3881 (7/17) Page 1 of 1 PSC Publishing Services (301) 445-0740 EF

No. 1 Road, Long Thanh Industrial Zone, Taman Village, Long Thanh District, Dong Nai Province, Vietnam Tel: 84-61-3514022

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2020 See PRA Statement below

510(k) Number (if known)

K190403

Device Name

Disposable Powder Free Nitrile Examination Glove, White Color, Tested For Use With Chemotherapy Drugs

Indications for Use (Describe)

A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner.

In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs:

Test Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Min)
1. Carmustine (BCNU), 3.3 mg/ml	22.8
2. Thiotepa 10.0 mg/ml	54.6
Cyclophosphamide (Cytoxan), 20.0 mg/ml	≥ 240
4. Cisplatin, 1.0 mg/ml	≥ 240
Doxorubicin Hydrochloride, 2.0 mg/ml	≥ 240
6. Fluorouracil, 50.0 mg/ml	≥ 240
Dacarbazine (DTIC), 10.0 mg/ml	≥ 240
8. Etoposide (Toposar). 20.0 mg/ml	≥ 240
9. Paclitaxel (Taxol), 6.0 mg/ml	≥ 240

The maximum testing time is 240 minutes. Please note that the following drug has an extremely low permeation time Carmustine (BCNU), 3.3 mg/ml

22.8 minutes

Thiotepa, 10.0 mg/ml 54.6 minutes

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

"DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW."

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRA Staff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

FORM FDA 3881 (7/17) Page 1 of 1

No. 1 Road, Long Thanh Industrial Zone, Taman Village, Long Thanh District, Dong Nai Province, Vietnam Tel: 84-61-3514022

510(k) SUMMARY

K190403

Premarket Notification [510(k)] Summary

1.0 Submitter:

Submitter's name : Ever Global (Vietnam) Enterprise Corporation

Submitter's address: No. 1 Road, Long Thanh Industrial Zone, Taman Village, Long Thanh District,

Dong Nai Province, Vietnam

 Phone number:
 84-61-3514022

 Fax number:
 84 -61-3514023

Name of contact person: Jerry Lin
Summary Preparation Date: September 6, 2019

2.0 Name of the Device

Proprietary/Trade name: Disposable Powder Free Nitrile Examination Glove, White Color, Tested For

Use With Chemotherapy Drugs, Disposable Powder Free Nitrile Examination

Glove, Blue Color, Tested For Use With Chemotherapy Drugs

Common Name: Nitrile Examination Gloves
Classification Name: Patient Examination Glove

Device Classification: Class I

Regulation Number: 21 CFR 880.6250

Product Code: LZA, LZC

3.0 Predicate device

Device Name: KS Medicare Powder Free Nitrile Examination Gloves, Non-sterile, Tested for

Use with Chemotherapy Drugs (Blue, Black)

Company name: Koon Seng Sdn. Bhd.

510(K) Number: K171171

4.0 Device Description:

Disposable Powder Free Nitrile Examination Glove, White/Blue Color, Tested For Use With Chemotherapy Drugs is a patient examination glove made from nitrile compound, non-sterile (as per 21 CFR 880.6250, Class I). The principle operation of the medical device to provide single use barrier protection for the wearer and the device meets all the requirement specifications for Barrier Protection, tensile properties as defined in ASTM D6319-10, Standard specification for Nitrile Examination Gloves.

5.0 Indication for use:

Disposable Powder Free Nitrile Examination Glove, White Color, Tested For Use With Chemotherapy Drugs

A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs:

Table 5.1 Nitrile white color test for 9 chemotherapy drugs

	Test Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Min.)	
1.	Carmustine (BCNU), 3.3 mg/ml	22.8	
2.	Thiotepa, 10.0 mg/ml	54.6	
3.	Cyclophosphamide (Cytoxan), 20.0 mg/ml	2 240	
4.	Cisplatin, 1.0 mg/ml	2 240	
5.	Doxorubicin Hydrochloride, 2.0 mg/ml	2 240	
6.	Fluorouracil, 50.0 mg/ml	2 240	
7.	Dacarbazine (DTIC), 10.0 mg/ml	2 240	
8.	Etoposide (Toposar), 20.0 mg/ml	2 240	
9.	Paclitaxel (Taxol), 6.0 mg/ml	2 240	
The maximum testing time is 240 minutes. Please note that the following drug has an			

The maximum testing time is 240 minutes. Please note that the following drug has an extremely low permeation time:

Carmustine (BCNU), 3.3 mg/ml	22.8 minutes	
Thiotepa, 10.0 mg/ml	54.6 minutes	ı

No. 1 Road, Long Thanh Industrial Zone, Taman Village, Long Thanh District, Dong Nai Province, Vietnam Tel:84-61-3514022

$\underline{Disposable Powder Free Nitrile Examination Glove, Blue Color, Tested For Use With Chemotherapy Drugs$

A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs:

Table 5.2 Nitrile blue color test for 9 chemotherapy drugs

	Test Chemotherapy Drug and Concentration	Minimum Breakthrough Detection Time (Min.)	
1.	Carmustine (BCNU), 3.3 mg/ml	6.2	
2.	Thiotepa, 10.0 mg/ml	38.8	
3.	Cyclophosphamide (Cytoxan), 20.0 mg/ml	2 240	
4.	Dacarbazine (DTIC), 10.0 mg/ml	2 240	
5.	Doxorubicin Hydrochloride, 2.0 mg/ml	2 240	
6.	Etoposide (Toposar), 20.0 mg/ml	2 240	
7.	Fluorouracil, 50.0 mg/ml	2 240	
8.	Paclitaxel (Taxol), 6.0 mg/ml	2 240	
9.	Cisplatin, 1.0 mg/ml	2 240	

The maximum testing time is 240 minutes. Please note that the following drug has an extremely low permeation time:

Carmustine (BCNU), 3.3 mg/ml	6.2 minutes	l
Thiotepa, 10.0 mg/ml	38.8 minutes	

6.0 Technological Characteristics:

Shown below is a technological comparison of the subject device(K190403) with the predicate device (K171171)

Device Characteristic	Predicate Device (K171171)	Proposed Device (K190403)	Comparison
Product name	KS Medicare Powder Free Nitrile Examination Gloves, Non-sterile, Tested for Use with Chemotherapy Drugs (Blue, Black)	Disposable Powder Free Nitrile Examination Glove, White/Blue Color, Tested For Use With Chemotherapy Drugs	N/A
510(K) No.	K171171	K171422	N/A
Product Owner	Koon Seng Sdn. Bhd.	Ever Global Enterprise Corporation	Different
Product Code	LZA, LZC	LZA, LZC	same
Regulation	21 CFR 880.6250	21 CFR 880.6250	same
Class	I	I	same
Intended Use	A powder-free patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs:	The Nitrile Powder Free patient examination glove is a non-sterile disposable device intended for medical purposes that is worn on the examiner's hands or finger to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs:	similar
Power free	Yes	Yes	same
Size	Small/ Medium/Large/X Large	Small/ Medium/Large/X Large	same
Single Use	YES	YES	same
Non-Sterile	YES	YES	same
Dimensions- Length	Complies with ASTM D6319-10 230 mm min.	Complies with ASTM D6319-10 230 mm min.	same

No. 1 Road, Long Thanh Industrial Zone, Taman Village, Long Thanh District, Dong Nai Province, Vietnam Tel:84-61-3514022

1011011011011						
Dimensions	Complies with ASTM D6319-10		Complies with A	Complies with ASTM D6319-10		
-Palm Width	Small	80±10	Small	80±10		
	Medium	95 ±10	Medium	95 ±10		
	Large	110±10	Large	110±10		
	X large	120 ±10	X large	120 ±10		
Dimensions	Complies with ASTM D6319-10		Complies with A	Complies with ASTM D6319-10		
-Thickness	-Thickness Palm - 0.07 mm min.		Palm - 0.05mi	Palm - 0.05mm min.		
	Finger - 0.08 mm min.		Finger - 0.05 mi	Finger - 0.05 mm min.		
			{Additional con	{Additional comparisons a device - K171104,		
			that have been	that have been approved by FDA is the same		
			intended use a	intended use and the thickness is lower than		
			Proposed Device	e)		

Intended Use	A consider for a matter to consider the con-	The Minister December Francisco Control Control	-11
intended Use	A powder-free patient examination glove is	The Nitrile Powder Free patient examination	similar
	a disposable device intended for medical	glove is a non-sterile disposable device	
	purposes that is worn on the examiner's	intended for medical purposes that is worn on	
	hand or finger to prevent contamination	the examiner's hands or finger to prevent	
	between patient and examiner.	contamination between patient and examiner.	
	In addition, these gloves were tested for use	In addition, these gloves were tested for use	
	with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for	with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for	
		Assessment of Medical gloves to Permeation	
	Assessment of Medical gloves to Permeation	S	
Power free	by Chemotherapy Drugs: Yes	by Chemotherapy Drugs: Yes	
Size			same
	Small/ Medium/Large/X Large	Small/ Medium/Large/X Large	same
Single Use	YES	YES	same
Non-Sterile	YES	YES	same
Dimensions- Length	Complies with ASTM D6319-10	Complies with ASTM D6319-10	same
•	230 mm min.	230 mm min.	
Dimensions -Palm Width	Complies with ASTM D6319-10	Complies with ASTM D6319-10	same
-raiiii vvidtii	Small 80±10	Small 80±10	
	Medium 95 ±10	Medium 95 ±10	
	Large 110±10	Large 110±10	
	X large 120 ±10	X large 120 ±10	
Dimensions	Complies with ASTM D6319-10	Complies with ASTM D6319-10	
-Thickness	Palm - 0.07 mm min.	Palm - 0.05mm min.	similar
	Finger - 0.08 mm min.	Finger - 0.05 mm min.	
		{Additional comparisons a device - K171104,	
		that have been approved by FDA is the same	
		intended use and the thickness is lower than	
		Proposed Device)	
Biocompatibility	Primary Skin Irritation in rabbits ISO 10993-10	Primary Skin Irritation in rabbits ISO 10993-10:	
	Dermal sensitization in the guinea pig ISO	Third Edition 2010-08-01	similar
	10993-10	Dermal sensitization in the guinea pig ISO 10993-	
	/	10: Third Edition 2010-08-01	
		Minimal Essential Media (MEM) Elution	
		ISO 10993-5	

7.0 Summary of Non-Clinical Performance Data:

Shown below are the results from the bench testing performed on the subject devices. The testing was performed to $\frac{1}{2}$

No. 1 Road, Long Thanh Industrial Zone, Taman Village, Long Thanh District, Dong Nai Province, Vietnam Tel:84-61-3514022

demonstrate the subject devices met the acceptance criteria for each of the non-clinical test shown below:

Disposable Powder Free Nitrile Examination Glove, White/Blue Color, Tested For Use With Chemotherapy Drugs are summarized with the following technological characteristics compared to ASTM or equivalent standard.

Table 7.1 Summary of the Non-Clinical Testing

Characteristics	Standard				
Dimension					
	Length	ngth 2230mm			
	Width	Small	<u>80 ±10 mm</u>		
		Medium	<u>95 ±10 mm</u>		
		Large	<u>110 ±10</u> mm		
		X large	<u>120 ±10</u> mm		
	Thickness	Finger tip	20.05mm		
		Palm	20.05mm		
Physical Properties	ASTM standard D 6319-	10(Reapproved 2015	5)		
	Tensile strength (Before	aging)		214MPa	
	Tensile strength (After ag	ging)		214MPa	
	Elongated rate (Before a	ging)		2500%	
	Elongated rate (After agi	ing)		2400%	
Freedom from pinholes	21 CFR 800.20			Passed Standard	
	ASTM standard D 6319-			Acceptance	
	Test method in accordan		· ' ' ' ' ' ' ' ' ' ' ' ' ' ' ' ' ' ' '	Criteria	
Powder Residual	ASTM standard D6319-3				
	Test method in accordance with D6124-06(R		,	< 2 mg/glove	
Biocompatibility	Primary Skin Irritation in rabbits		Passes		
	ISO 10993-10: Third Edition 2010-08-01		Under the conditions of the study, the subject		
	Dermal sensitization in the guinea pig ISO 10993-10: Third Edition 2010-08-01		device is not a primary skin irritant.		
			Passes		
			Under the conditions of the study, the subject		
			device is not a primary skin sensitizer.		
	Minimal Essential Media ISO 10993-5	(IVIEIVI) EIUTION	Passes		
	130 10332-3		Under the conditions of the study, the subject device is not cytotoxic.		

- Dimension per ASTM D6319-10 (Reapproved 2015)
- ★ Water leak test on pinhole per ASTM D6319-10(Reapproved 2015) and per 21 CFR 800.20.
- Biocompatibility test per ISO 10993-10: Third Edition 2010-08-01 and ISO 10993-5
- ☆ Assessment of Resistance To Permeation By Chemotherapy Drugs per ASTM D6978-05(R 2013)

8.0 Summary of Clinical Test:

Clinical data was not needed for this device.

9.0 Conclusion:

Based on the nonclinical tests data, it can be concluded that the Disposable Powder Free Nitrile Examination Glove, White/Blue Color, Tested For Use With Chemotherapy Drugs is as safe, as effective, and performs as well as or better than the predicate device K171171, KS Medicare Powder Free Nitrile Examination Gloves, Non-sterile, Tested for Use with Chemotherapy Drugs (Blue, Black), by Koon Seng Sdn. Bhd.

CE CERTIFICATION



CE Registration Certificate

This is to certify that, in accordance with the Medical Device Directive 93/42/EEC, Emergo Europe agrees to perform all duties and responsibilities as the Authorized Representative for

as stipulated and demanded by the aforementioned Directive. The Dutch Competent Authorities have received Medical Device Registrations on the following date:

> 26 May 2009 See attached product listing

Emergo Europe Registration Number: NL/CA01/601529

The Manufacturer has provided Emergo Europe with the appropriate Declaration (Manufacturer has provided Emergo Europe with the applicable requirements of Directive 93/42/EEC.

June 2009

Rene van de Zande President & CEO Emergo Europe

CERTAND STANDSTAND

Emergo Europe





Technical File	Doc. No.	PM-TCF-DOC-PPE		
recunical rile	Rev. No.	1B		
EU 2016/425	Rev. Date	2018.11.09		
EU 2016/425	Page	1/1		

EU Declaration of Conformity

We, **PRECIOUS MOUNTAIN ENT. CORP.**, hereby declares that the following Personal Protective Equipment (PPE).

1. Model & Description

■ GNP (Disposable Nitrile Gloves, Powdered)

■ GNF (Disposable Nitrile Gloves, Powder-free)

2. Glove sizes

XS, S, M, L, XL or 6, 7, 8, 9, and 10

3. Manufacturer

PRECIOUS MOUNTAIN ENT. CORP.

2F No.68, Sec 1., Neihu rd., Neihu Dist., Taipei City 114, Taiwan

P. 886-2-6606-6688 P. 886-2-6606-1818

4. Classification

Category III-Complex

- 5. It is in conformance with the Regulation (EU) No 2016/425 of the European Parliament and of the Council and national standards in some cases, in conformance with the harmonized standards conversion.
- 6. The notified body SATRA (Number: 2777) performed the EU type-examination certificate number is_____ (Module B).
- 7. The object of the declaration described in point 4 is in conformity with the relevant Union harmonization legislation: EN420:2003+A1:2009, EN374-1:2016, EN ISO 374-5:2016.
- 8. Europe Standard
- ✓ EN ISO 374-1:2016,
- ✓ EN 374-2: 2014
- ✓ EN 16523-1:2015
- ✓ EN 420:2003+A1:2009
- ✓ EN ISO 374-4:2013
- ✓ EN ISO 374-5:2016,

EN ISO 374-1:2016/Type B



(B) (i

EN ISO 374-5: 2016

(E

For and on behalf of

PRECIOUS MOUNTAIN ENT. CORP.

Tammy Phone / Quality Manager

Date: Nov. 09, 2018

TEST REPORT



Signature Not Verified For Question Please Contact with SGS www.sgs.com.tw

Test Report No.: CT/2018/70053C

Date: 2018/09/18

Page: 1 of 4

PRECIOUS MOUNTAIN ENT. CORP

2F., NO. 68, SEC. 1, NEIHU RD., NEIHU DIST., TAIPEI, TAIWAN

The following sample(s) was/were submitted and identified by/on behalf of the applicant as:

Sample Submitted By

: PRECIOUS MOUNTAIN ENT. CORP

Sample Description

: DISPOSABLE NITRILE GLOVES (BLUE)

Style/Item No.

: LOT#HL12050401

Sample Material

: NBR(N300ITRILE)

Sample description(Item No.) : DISPOSABLE NITRILE GLOVES (BLACK/WHITE/VIOLET)

of the same material

Sample Receiving Date

: 2018/07/03

Testing Period

2018/07/03 TO 2018/08/02 AND 2018/07/19 TO 2018/08/02 AND 2018/08/16 TO 2018/08/27

Test Requested

: As specified by client, the sample(s) was/were tested with reference to Resolution (EC) No 1935/2004, Council of Europe Resolution ResAP(2004)4. Please refer to result table for

testing item(s).

Test Result(s)

Signed for and on beh SĞS TAIWAN LTD

Chemical Laboratory - Taipei

: Please refer to following pages.

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <a href="http://www.sgs.com/en/Terms-and-Conditions/ferms-

^{*} This report is added testing and combined with CT/2018/70053 *



Test Report No.: CT/2018/70053C

Date: 2018/09/18

Page: 2 of 4

PRECIOUS MOUNTAIN ENT. CORP

2F., NO. 68, SEC. 1, NEIHU RD., NEIHU DIST., TAIPEI, TAIWAN

Test Result(s)

PART NAME No.1

: BLUE RUBBER GLOVE

PASS

Test Item(s)	Unit	Method	MDL	Result No.1	Limit
Overall migration (3% Acetic Acid)	mg/dm²		3.0	7.83	10
Overall migration (10% Ethanol)	mg/dm²	With reference to EN 1186 (2002). (70°C, 2 hrs)	3.0	9.50	10
Overall migration (Rectified Olive Oil)	mg/dm²	(70 C, 2 hrs)	3.0	n.d.	10
Nitrosamines					
N-nitrosodimethylamine (CAS No.: 62-75-9)	mg/kg		0.01	n.d.	0.01
N-nitrosodiethylamine (CAS No.: 10595-95-6)	mg/kg		0.01	n.d.	0.01
N-nitrosomethylethylamine (CAS No.: 55-18-5)	mg/kg		0.01	n.d.	0.01
N-nitrosodipropylamine (CAS No.: 621-64-7)	mg/kg		0.01	n.d.	0.01
N-nitrosomorpholine (CAS No.: 59-89-2)	mg/kg		0.01	n.d.	0.01
N-nitrosopyrrolidine (CAS No.: 930-55-2)	mg/kg	With reference to EN 14350 (2004).	0.01	n.d.	0.01
N-nitrosopiperidine (CAS No.: 100-75-4)	mg/kg	Analysis was performed by GC/MS.	0.01	n.d.	0.01
N-nitrosodibutylamine (CAS No.: 924-16-3)	mg/kg		0.01	n.d.	0.01
N-nitrosodiphenylamine (CAS No.: 86-30-6)	mg/kg		0.01	n.d.	0.01
Methylphenylnitrosamine (CAS No.: 614-00-6)	mg/kg		0.01	n.d.	0.01
Ethylphenylnitrosamine (CAS No. 612-64-6)	mg/kg		0.01	n.d.	0.01
N-nitrosodibenzylamine (CAS No.: 5336-53-8)	mg/kg		0.01	n.d.	0.01
Nitrosatable					
N-nitrosodimethylamine	mg/kg		0.1	n.d.	0.1
N-nitrosodiethylamine	mg/kg		0.1	n.d.	0.1
N-nitrosomethylethylamine	mg/kg		0.1	n.d.	0.1
N-nitrosodipropylamine	mg/kg		0.1	n.d.	0.1
N-nitrosomorpholine	mg/kg		0.1	n.d.	0.1
N-nitrosopyrrolidine	mg/kg	With reference to EN 14350 (2004).	0.1	n.d.	0.1
N-nitrosopiperidine	mg/kg	Analysis was performed by GC/MS.	0.1	n.d.	0.1
N-nitrosodibutylamine	mg/kg		0.1	n.d.	0.1
N-nitrosodiphenylamine	mg/kg		0.1	n.d.	0.1
Methylphenylnitrosamine	mg/kg		0.1	n.d.	0.1
Ethylphenylnitrosamine	mg/kg		0.1	n.d.	0.1
N-nitrosodibenzylamine	mg/kg		0.1	n.d.	0.1

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <a href="http://www.sus.com/en/Terms-and-Conditions/Terms-



Test Report No.: CT/2018/70053C

Date: 2018/09/18

Page: 3 of 4

PRECIOUS MOUNTAIN ENT. CORP

2F., NO. 68, SEC. 1, NEIHU RD., NEIHU DIST., TAIPEI, TAIWAN

THE RESERVE OF	norman a	THE REAL PROPERTY.	monenn

Test Item(s)	Unit	Method	MDL	Result No.1	Limit
Specific migration of Primary Aromatic Amine (PAA)	mg/kg	With reference to EN 13130-1 (2004) and BfR recommendation XXI. Analysis was performed by UV-vis Spectrophotometer. (3% acetic acid, 70°C, 2 hrs)	0.002	n.d.	0.01 (▲1)

Remark:

- (A1): Limit is according to Council of Europe Resolution AP (2004) 4.

Note:

- 1. mg/kg = ppm; mg/dm² = milligram per square decimeter
- 2. MDL = Method Detection Limit
- 3. n.d. = Not Detected = below MDL
- 4. °C = degree Celsius
- 5. Test condition & simulant were specified by client.
- 6. This report supersedes the previous document bearing the test report number CT/2018/70053A.

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <a href="http://www.sgs.com/en/Terms-and-Conditions-forms-and-Conditions-and-Conditions-and-Conditions-forms-and-Conditio



Test Report No.: CT/2018/70053C

Date: 2018/09/18

Page: 4 of 4

PRECIOUS MOUNTAIN ENT. CORP 2F., NO. 68, SEC. 1, NEIHU RD., NEIHU DIST., TAIPEI, TAIWAN

* The tested sample / part is marked by an arrow if it's shown on the photo. *

CT/2018/70053



** End of Report **

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-en-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of client's instruction, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced, except in full, without prior written approval of the Company, Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the sample(s) tested.

ISO 13485:2016 & 9001:2015





Certificate

No. Q6 098755 0002 Rev. 00

Holder of Certificate: Ever Global (Vietnam) Ent. Corp -

Giang Dien Factory

No 9 Street, Giang Dien Industrial Park

Trang Bom District Dong Nai Province

VIETNAM

Ever Global (Vietnam) Ent. Corp - Giang Dien Factory Facility(ies):

No 9 Street, Giang Dien Industrial Park, Trang Bom District, Dong

Nai Province, VIETNAM

Certification Mark:



Scope of Certificate: Manufacturing and sales of powdered and

powder free synthetic vinyl gloves and powder free nitrile gloves for medical use

Applied Standard(s): EN ISO 13485:2016

Medical devices - Quality management systems -

Requirements for regulatory purposes

(ISO 13485:2016) DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system (excluding subclause 7.3). which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.: VIE201802

Valid from: 2018-07-26 Valid until: 2020-05-14

2018-07-26 Date,

1. Pumil



CERTIFICATE

The Certification Body of TÜV SÜD Management Service GmbH

certifies that

Ever Global (Vietnam) Ent. Corp - Giang Dien Factory No 9 Street, Giang Dien Industrial Park, Trang Bom District, 70000 Dong Nai Province Vietnam

> has established and applies a Quality Management System for

Manufacturing and sales of powdered and powder free synthetic vinyl gloves and powder free nitrile gloves for medical use.

An audit was performed, Report No. 707077466.

Proof has been furnished that the requirements according to

ISO 9001:2015

are fulfilled.

The certificate is valid from 2018-08-10 until 2020-07-23.

Certificate Registration No.: 12 100 54329 TMS.

Product Compliance Management Munich, 2018-08-13













SUPÉJIEUT® SIMPLY SUPERIOR

Nitrile Examination Gloves

• Powder Free • Ambidextrous • Single Use Only • Non-sterile

10 Boxes 100 Per Box







☑ Blue

☐ Black

☐ White

EVER GLOBAL (VIETNAM) ENTERPRISE CORP.

Made in Vietnam







Nitrile Examination Gloves

Nitrile Butadiene Rubber Powder free, Non-sterile, Ambidextrous. Free from natural rubber latex.







Storage: Protect from heat, humidity, strong light and ozone.



LOT 20200213



022020



022024

Made in Vietnam

EVER GLOBAL (VIETNAM) ENTERPRISE CORPORATION

Address : Route 1, Long Thanh Industrial Zone, Long Thanh District,

Dong Nai Province, 76211, Vietnam

E-mail: info@egvnco.com


TEL: 0251-8966-677

License No.: 190000024/PCBA-DN 180000010/PCBSX-DN

Item No. : NBR027



Select Your Glove Size







⊗ (€ 🛍 www.pmgloves.com





Nitrile

Powder Free
Nitrile Examination Gloves S

